

Snap-On wire markers

Tyco Electronics snap-on markers are convenient and easy to use. These markers are supplied on a 'wand' applicator that simplifies and reduces installation time, and saves on the cost of marking.

Snap-on wire markers are designed to stand the test of time, preventing costly re-work before the life of the equipment ends. This is why Tyco Electronics snap-on markers are found today in power stations, oil refineries, offshore platforms, and other major projects around the world.


Features and benefits

- Snap-on markers are the perfect way to permanently and rapidly mark cables, wires, leads, and cores after termination and installation
- Easy to alter identification on cables that are already marked
- Secure grip for positive fixing to cable
- The high elasticity of the marker enables all types of installation and conduction material to be permanently marked without damage
- Applicator has a patented slit to prevent over-expansion of the marker during application
- Made from polyoxymethylene which guarantees strength, dimensional stability, low moisture absorption, and resistance to chemicals
- Resistant against oils, cleaning agents, fuels, and alkaline solutions
- Zero-halogen material
- Ten sizes for cable diameters 0.04" to 0.75" (1.0mm to 19.0mm)
- Color-coded option



| TDS-094 | |
|----------------|-----------------|
| | |
| 40°C to +106°C | -40°F to +223°F |
| | |

Part numbering system









Snap-On wire markers

Ordering information

| Size Code | Cable D mm | Diameter inches | Color of Applicator Wand | Packed Box |
|--------------|---------------|-----------------|--------------------------|---|
| STD 01 | 1.00 - 1.40 | 0.040 - 0.055 | Brown | 10 mounting wands with 30 identical markers |
| STD 02 | 1.40 - 1.80 | 0.055 - 0.071 | Gray | 10 mounting wands with 30 identical markers |
| STD 03 | 1.90 - 2.65 | 0.075 - 0.104 | Green | 10 mounting wands with 30 identical markers |
| STD 05 | 2.30 - 3.00 | 0.091-0.118 | Natural | 10 mounting wands with 30 identical markers |
| STD 06 | 2.60 - 3.50 | 0.102 - 0.138 | Red | 10 mounting wands with 30 identical markers |
| STD 09 | 3.25 - 4.50 | 0.128 - 0.177 | Blue | 10 mounting wands with 30 identical markers |
| STD12 | 4.50 - 6.00 | 0.177-0.236 | Yellow | 10 mounting wands with 30 identical markers |
| STD 15 | 5.80 - 8.50 | 0.228 - 0.335 | Orange* | 50 identical markers; applicator available on request |
| STD 17 | 8.50 - 11.50 | 0.335 - 0.452 | Orange* | 50 identical markers; applicator available on request |
| STD 21 | 11.00 - 15.50 | 0.433 - 0.610 | Gray* | 50 identical markers; applicator available on request |
| STD 24 | 15.00 - 19.00 | 0.590 - 0.750 | Gray* | 50 identical markers; applicator available on request |
| | | | | |

Color Code (Available for All Legends)

Y Yellow

W White

Leaend code

| Legena cod | ie | | | | | |
|------------|------|--------|-------|--------|--------|--|
| Legend | Code | Legend | Code | Legend | Code | |
| 0 | 0 | J | J | / | Slash | |
| 1 | 1 | K | K | Earth | Earth | |
| 2 | 2 | L | L | Blank | Blank | |
| 3 | 3 | M | M | : | Colon | |
| 4 | 4 | N | N | = | Equals | |
| 5 | 5 | 0 | 0 | | Dot | |
| 6 | 6 | Р | Р | | | |
| 7 | 7 | Q | Q | | | |
| 8 | 8 | R | R | | | |
| 9 | 9 | S | S | | | |
| А | А | T | T | | | |
| В | В | U | U | | | |
| С | С | V | V | | | |
| D | D | W | W | | | |
| Е | E | X | X | | | |
| F | F | Υ | Υ | | | |
| G | G | Z | Z | | | |
| Н | Н | + | Plus | | | |
| 1 | 1 | - | Minus | | | |
| | | | | | | |

Applicators

| | P.P. Commission | | | |
|--|---------------------|-------------|---------------------|------------------------------------|
| | Order Code | Marker Size | Color of Applicator | Number of Applications Per Pack |
| | STD 15 - Applicator | 15 | Orange | 5 |
| | STD 17 - Applicator | 17 | Orange | 5 |
| | STD21-Applicator | 21 | Gray | 2 |
| | STD24- Applicator | 24 | Grav | 2 |

83